




Controlling Cost And Quality Of

VITREORETINAL

Surgery*

* Maximizing pharmacoeconomic value (best health outcome for the resources invested)



HASA OPTIX
Belgium

Build your own vitreoretinal set



Tailor-made: Build your own instrument set from a portfolio of 1200 instruments which allows surgeons to prepare ahead of time for success

Simplify: Reduce number of SKU's which contributes a better work environment for all stakeholders

Constant quality: HASA Optix offers a quality equivalent to reusable devices allowing surgeons to use them constantly

Environmentally friendly: Instrument sets are presented in a paper pouche avoiding heavy plastic inserts (only one pack and one sterilization)


Circular economy: We are building a recycling program to give a second life to the stainless steel our instruments are made of

Cost controlling: We are offering intergrated solutions to allow the calculation of a fixed cost per surgical act


Include your vitreoretinal instrument sets in our vitreoretinal procedure pack for better OR organization and cost controlling



Anatomy Forceps



Anatomy Speculum



Anatomy Micro instruments



SCISSORS



Curved -
Horizontal blades




Straight -
Horizontal blades



Vertical -
Vertical blades

FORCEPS



End Grasping -
2.3mm long jaws



ILM -
1.25mm long jaws



Serrated-
1.25 long jaws



ICL -
Atraumatic jaws with holes



Capsulorhexis -
Angled grasping platform

PUNCHES



Punch -
Trabeculectomy
punch

19 Gauges

20 Gauges

23 Gauges

25 Gauges

27 Gauges

HASA's surgeons most wanted


Speculums



655035




Swivel



555025



655041



555023



655040

655035 Lieberman K wire adjustable blades

555025 Closed wire adjustable blades

655041 Solid round adjustable blades

555023 Lieberman V wire adult adjustable

655040 Open wire round adjustable blades

Indenter



554026



654024

554026 Schocket Scleral Depressor

654024 Indenter


Forceps



651080



551001



551175



551187

Scissors




653034



653055




553014



553016

Needle holders



552021



552022



552029



552007

Retina



651080 Paufigue Forceps

551001 Bonn Forceps 0.12 – Flat Handle 85 mm

551175 Tying Forceps - Curved jaws, flat handle 85 mm

551187 Colibri forceps 0.12 – Serrated handle

653034 Westcott Scissors Blunt Curved

653045 Westcott Scissors Blunt Straight

653055 Westcott Scissors Sharp Curved

653056 Westcott Scissors Sharp Straight

553014 Vannas curved sharp blades

553016 Vannas straight sharp blades

552021 9 mm curved jaws Round knurled handle

552022 10 mm curved jaws Round knurled handle

552029 9 mm curved jaws Flat serrated handle

552007 9 mm straight jaws Flat serrated handle

Retina Forceps

556038 23/25G Eckardt end-gripping

556026 23/25G ILM Forceps

556030 23/25G Serrated forceps 1.25 mm long jaws

Retina Scissors

556006 23/25G Horizontal curved scissors 1.25 mm blades

556003 23/25G Horizontal straight scissors 1.25 mm blades

556010 23/25G Vertical scissors 80°

556038

556026

556030

556006

556003

556010

Cataract

Cornea

Refractive

Glaucoma

Oculoplastic

Injection


Vitreoretinal

Anticipate, simplify and professionalize are our drivers.
Join us to help rethink ophthalmic surgery

Contact us to find your nearest distributor



 +32.2.524.63.88

 Info@hasaoptix.com

 www.hasaoptix.com

 Follow us

Avenue Louise 523, 1050 Brussels
Belgium

